

PROLAQ WS starter set

Item no. 837236

Detailed information on page 165-167



Scope of delivery:

- 1 x 837229 PROLAQ WS sedimentation tank
- 1 x 837230 PROLAQ WS
- 1 x 837235 PROLAQ WS Rinse

PROLAQ WS sedimentation tank

Item no. 837229

FOR SEDIMENTATION OF PROLAQ WS.

The sedimentation insert can be completely disassembled and thus enables easy cleaning of the system.

Scope of delivery:


- 1 x Stainless steel tank (30 l)
- 1 x sedimentation insert (5 pcs.)



PROLAQ WS

Item no. 837230

For the removal of solvent-based and water-based paints and varnishes



PROLAQ WS Rinse

Item no. 837235

For the rinsing of painting tools after cleaning with PROLAQ WS.

After cleaning with PROLAQ WS cleaners, rinse the painting tools manually in an immersion bath with PROLAQ WS Rinse. For machine use, add the rinsing liquid to the PROLAQ cleaning system.



CLEAN

BLOW DRY

RINSE

BLOW DRY



PROLAQ WS

PROLAQ WS Rinse